

DISCOVER for Crime-lite[®] AUTO

Maximise the search, detection, and imaging of evidence at the crime scene <u>and</u> in the laboratory

Combining the handheld Crime-lite AUTO with a modular laboratory dock, DISCOVER is a versatile two-part workstation for every stage of the forensic evidence workflow.



Un-docked, the **Crime-lite AUTO** is a complete CSI search and capture tool for the multi-spectral detection and imaging of evidence at the crime scene.

Docked, the **DISCOVER** workstation adds specialist imaging and illumination modules for the photography of fingermarks. Dedicated software provides casework management and report builder functions.

DISCOVER takes the intuitive operation of the Crime-lite AUTO and applies it to laboratory fingerprint imaging

Versatile Fingerprint Imaging System

-----Applications------

- Forensic Laboratory

Visualise and Capture Latent, Patent, & Processed Fingermarks:

Multi-Spectral Fluorescence Imaging				
Ultraviolet		Infrared		
Contrast	Coaxial		ial	Oblique
Cylindrical Surfaces			Polarised Light	

-Hardware---

Laboratory Docking System

Core System

- Crime-lite AUTO quick-dock port
- XYZ positioning
- Mouse/Keyboard/Monitor interface

Illumination & Imaging Modules

- Cylindrical Surface Unwrapper
- Coaxial Illumination
- Polarised Illumination
- Diffuse Light Plate
- Seamless transition between Crime Scene & Laboratory

-----Features-----

- Perform visible, fluorescent, contrast, UV and IR examinations
- Advanced software functionality



foster+freeman



Crime-lite® AUTO

The most complete and compact forensics digital camera solution in existence

Combining the latest forensic imaging technology with high-intensity multi-spectral illumination, the new Crime-lite AUTO is a complete solution to the search, detection, and capture of evidence.

A compact handheld device, the Crime-lite AUTO includes a 20MP camera, full-spectrum LED illumination, and provides intelligent filter selection to enable the examiner to quickly and easily reveal a wide range of evidence types.

Unlock the Crime-lite AUTO's full potential through the addition of a DISCOVER modular laboratory docking system (see over for details).


Forensic imaging technology meets powerful multi-spectral illumination, in this complete solution to the search, detection, and capture of evidence

Crime Scene Examination Without Limitations

-----Applications----

- Crime Scene - Forensic Laboratory

Search, Detect and Image:



-----Hardware-----

- 20MP Camera
- UV-Vis-IR Imaging
- 55x illumination wavebands
- Automated filter selection
- Auto-Focus capture
- Touchscreen display
- Intuitive Interface
- -----Features------
- A complete, compact forensic search tool
- Real-time visible, fluorescent, and infrared examinations
- Search Beyond the Visible to detect every trace of evidence
- Optional laboratory docking system



foster+freeman UK info@fosterfreeman.com +44 (0)1386 768 050 foster+freeman USA usoffice@fosterfreeman.com 888 445 5048

foster+freeman EU de.sales@fosterfreeman.com +49 (0)211 158 422 57 foster+freeman NL nl.sales@fosterfreeman.com +31 (0)6 1114 3858 foster+freeman