

VSC[®] EXTERNAL MICROSCOPES

Combine the document examination capabilities of the VSC[®]8000 with one of our selected microscopes to achieve the ideal interaction of multi-spectral imaging and advanced image processing with crystal clear imaging of fine details at high levels of magnification.



PVM Portable Video Microscope

VSC/PVM

This compact 3x optical zoom video microscope with integrated LED illumination, provides a premium, portable solution to the examination of taggants.

Supplied with a 5MP USB CMOS camera and stand-alone software, the PVM's clear base allows easy positioning on documents and for the use of VSC[®]8000 oblique lighting.

3: 1

Zoom

249x

Magnification*

Vis LED

Illumination



DVM Digital Video Microscope

V6/DVM/3

Comprising a high specification microscope, 14 Bit CCD colour Firewire camera, two high intensity LED light sources (White and UV), darkfield ringlight and a transmitted white light XY stage, the DVM can be fully integrated with a VSC[®] system or used as a PC driven stand alone instrument complete with dedicated software.

10: 1

Zoom

10200x

Magnification*

UV-Vis LED

Illumination