



Versatile Fingerprint Imaging System

-----Applications-----

- Forensic Laboratory

Visualise and Capture Latent, Patent, & Processed Fingermarks:

- Multi-Spectral Fluorescence Imaging
- Ultraviolet Infrared
- Contrast Coaxial Oblique
- Cylindrical Surfaces Polarised Light

-----Hardware-----

Laboratory Docking System

Core System

- Crime-lite AUTO quick-dock port
- XYZ positioning
- Mouse/Keyboard/Monitor interface

Illumination & Imaging Modules

- Cylindrical Surface Unwrapper
- Coaxial Illumination
- Polarised Illumination
- Diffuse Light Plate

-----Features-----

- Seamless transition between Crime Scene & Laboratory
- Perform visible, fluorescent, contrast, UV and IR examinations
- Advanced software functionality

DISCOVER for Crime-lite® AUTO

Maximise the search, detection, and imaging of evidence at the crime scene and in the laboratory

Combining the handheld Crime-lite AUTO with a modular laboratory dock, DISCOVER is a versatile two-part workstation for every stage of the forensic evidence workflow.



Un-docked, the Crime-lite AUTO is a complete CSI search and capture tool for the multi-spectral detection and imaging of evidence at the crime scene.



Docked, the DISCOVER workstation adds specialist imaging and illumination modules for the photography of fingermarks. Dedicated software provides casework management and report builder functions.

DISCOVER takes the intuitive operation of the Crime-lite AUTO and applies it to laboratory fingerprint imaging





Crime Scene Examination Without Limitations

-----Applications-----

- Crime Scene
- Forensic Laboratory

Search, Detect and Image:

Gunshot Residues	Blood	Fingerprints
Accelerants	DNA	Body Fluids
Bruising	Bone & Tooth	Shoeprints
Hair & Fibres	and more...	

-----Hardware-----

- 20MP Camera
- UV-Vis-IR Imaging
- 55x illumination wavebands
- Automated filter selection
- Auto-Focus capture
- Touchscreen display
- Intuitive Interface

-----Features-----

- A complete, compact forensic search tool
- Real-time visible, fluorescent, and infrared examinations
- Search Beyond the Visible to detect every trace of evidence
- Optional laboratory docking system

Crime-lite® AUTO

The most complete and compact forensics digital camera solution in existence

Combining the latest forensic imaging technology with high-intensity multi-spectral illumination, the new Crime-lite AUTO is a complete solution to the search, detection, and capture of evidence.

A compact handheld device, the Crime-lite AUTO includes a 20MP camera, full-spectrum LED illumination, and provides intelligent filter selection to enable the examiner to quickly and easily reveal a wide range of evidence types.

Unlock the Crime-lite AUTO's full potential through the addition of a DISCOVER modular laboratory docking system (see over for details).

Forensic imaging technology meets powerful multi-spectral illumination, in this complete solution to the search, detection, and capture of evidence

